Version

1

CARDS MANUAL

Neuropsychological and Clinical Examination for Clinical and Population Use

Cognitive Assessment & Reference Diagnosis System (CARDS)

NEUROPSYCHOLOGICAL AND CLINICAL EXAMINATION

Cognitive Assessment & Reference Diagnosis System

© Guk-Hee Suh, M.D., Ph.D.
Department of Psychiatry
Hallym University Medical Center
Hangang Sacred Heart Hospital
95-200 Yungdungpo-Dong
Seoul, 150-030
Korea
Phone 82.2.2639.5483 • Fax 82.2.2677.9095
E-mail : suhgh@chollian.net

All Rights Reserved. No part of this publication may be reproduced, stored in a retrieval system or transmitted, in any form or by any means, electronic, mechanical, photocopying, recording, scanning or otherwise, except under the terms of the Copyright, Designs and Patents Act in Korea or under the terms of a licence issued by the Copyright Holder, without the permission in writing of the Copyright holder, Dr. Guk-Hee Suh, MD, PhD. Department of Psychiatry, Hallym University Medical Center, Hangang Sacred Heart Hospital, 95-200 Yungdungpo-Dong, Yungdungpo-Ku, Seoul, Korea (Zipcode: 150-030). E-mail : suhgh@chollian.net

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2003. 8. 6

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CARDS?

A new computerized neuropsychological interview and clinical examination to incorporate in a single instrument all the information required to make an accurate clinical diagnosis of dementia, its subtype, depression and delirium.

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2. Interview with subject
3. Interview with informants
4. interviewer's judgment
5. clinician's examination
6. laboratory tests and brain imaging

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Table of Contents

PREFACE	I
CARDS?	III
WHY CARDS?	1
WHY NOW?	2
WHAT'S THE PURPOSES?	3
BRIEF OUTLINE.....	3
VALIDITY & RELIABILITY OF THE CARDS	5
GENERAL INSTRUCTION FOR INTERVIEWER.....	7
HOW TO DEAL WITH PROBLEMS RELATED TO AGEING AND PHYSICAL ILLNESSES?	9
INTERVIEWER'S JUDGMENT	9
QUALIFYING INTERVIEWER	10
CONTENTS.....	11
NEUROPSYCHOLOGICAL TEST BATTERY: CARDS-COG.....	17
<i>Orientation</i>	17
<i>Word List Memory-Immediate Recall</i>	20
<i>Object Naming Test</i>	21
<i>Remote Memory</i>	23
<i>Recent Memory</i>	24
<i>Word List Memory-Delayed Recall</i>	25
<i>Word List Memory -Recognition</i>	26
<i>Object List Memory-Delayed Recall</i>	27
<i>Object List Memory-Recognition</i>	28
<i>Identification Test</i>	29
<i>Combined Arithmetic Test</i>	30
<i>Serial 7 subtraction & Backward Citation of Weekdays</i>	30
<i>Ideational Praxis</i>	31
<i>Verbal Fluency</i>	32
<i>Word Definition</i>	34
<i>Constructional Praxis</i>	34
<i>Weather Today</i>	36
<i>Writing</i>	36
<i>Ideomotor Praxis</i>	37
<i>Abstract Thinking</i>	38
<i>Ideational Fluency</i>	39
<i>Visual Reasoning</i>	41
<i>Trail Making Test A</i>	42
<i>Face Identification</i>	43
<i>Object Constancy</i>	43
<i>Job Identification</i>	44
<i>Sign Identification</i>	45

<i>Judgment</i>	46
<i>Clock Face Reading</i>	47
<i>What time is it now?</i>	47
INTERVIEW WITH SUBJECT	48
INTERVIEWER'S JUDGMENT	54
<i>Spoken Language Ability</i>	54
<i>Listening Comprehension</i>	55
<i>Paraphasia, Word Finding Difficulty</i>	56
<i>Retardation</i>	57
<i>Depression</i>	58
<i>Flat affect</i>	58
<i>Perseveration</i>	58
<i>Incoherence/ Irrelevance</i>	58
INTERVIEW WITH INFORMANTS	59
CLINICIAN'S EXAMINATION	65
LABORATORY TESTS AND BRAIN IMAGING	67
OPERATING MANUAL	68
<i>First Screen</i>	68
<i>Beginning a new case</i>	70
<i>Interrupting an interview</i>	72
<i>Resuming an interview</i>	73
<i>Reports for a single subject</i>	75
<i>How to print Reports</i>	80
<i>How to see Reports later</i>	80
<i>How the interview ends</i>	82
<i>Correcting earlier answers</i>	82
<i>Correcting client's information</i>	82
<i>Comparing repeated measurements</i>	82
ITEMS FOR CLINICIAN'S JUDGMENT	84
INDEX	86

Why CARDS?

Dementia is a progressive syndrome that results in significant impairment to cognitive function and has a great social impact.

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Why now?

Recent therapeutic advances to delay the progression of dementia have renewed calls to improve the early detection and management of individuals with dementia.

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- **Interview with informants**

Longitudinal comparisons of various aspects of cognition such as memory, emotion, behaviour, and personality

Past history

Functional level of subject

- **Interviewer's judgment**

Emotion, behavior, personality and motion

Table 1. CARDS-cog subscales and sections

Subscale	Sections
Memory test	Word List Memory – Immediate Recall Word List Memory – Delayed Recall Word List Memory – Recognition Object List Memory – Delayed Recall Object List Memory – Recognition Recent Memory Remote Memory Orientation Today's Weather Current Time
Agnosia test	Object Naming Test Identification Test Object Constancy Test
Aphasia test	Object Naming Test* Verbal Fluency Word Definition Expressive Language Comprehension of Spoken Language Word-finding difficulty
Apraxia test	Ideational Praxis Ideomotor Praxis Constructional Praxis Writing
Dysexecution test	Abstract Thinking Ideational Fluency Visual Reasoning Judgment
Attention/calculation test	Trail Making Test A Serial 7 Subtraction Backward Citation of Weekdays Combined Arithmetic Test

*Both Agnosia test and Aphasia test are including Object naming test, so that its score is counted twice.

- **clinician's examination**

Physical examination, neurological examination and mental state examination

- **laboratory tests and brain imaging**

Physical condition and possible causes of cognitive decline

- **Reference diagnosis**

All the possible diagnoses based on the information during whole CARDS interview

- **DSM-IV diagnosis**

Diagnosis derived from the diagnostic algorithm in the CARDS according to DSM-IV diagnostic criteria

Validity & Reliability of the CARDS

CARDS is a valid and reliable instrument not only for clinical diagnosis, but also for the evaluation of impaired cognitive domains.

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CHAPTER 1. INTRODUCTION

Table 2. Sensitivity and Specificity of the diagnosis of dementia, Alzheimer's disease, vascular dementia and depression according to the CARDS when compared to the clinician's diagnosis as a gold standard

CARDS Diagnoses	2 X 2 Table			Sensitivity	Specificity	Correct Classification
		G.S. +	G.S. -			
	CARDS+	A	B			
	CARDS-	C	D			
Dementia		48	6	0.96	0.87	0.92
		2	40			
Alzheimer's disease		35	2	0.97	0.97	0.97
		1	58			
Vascular dementia		13	4	0.95	0.95	0.95
		1	78			
Depression		11	0	0.92	1.00	0.99
		1	84			

G.S. ; Gold standard, equivalent to clinician's diagnosis

Sensitivity=A/(A+C), Specificity=D/(B+D), Correct Classification=(A+D)/(A+B+C+D)

(adopted from Suh & Lee (2003), *Psychogeriatrics*)

Table 3. Mean scores of CARDS-cog total and subscales in GDS stage 3 group and GDS stage 4 group

Category	GDS Stage 3 (n=21)		GDS Stage 4 (n =37)		Difference	
	Mean	S.D.	Mean	S.D.	F value	P value
CARDS-cog total *	127.4	17.1	81.3	39.3	32.3	0.000
- Amnesia subscale*	49.4	7.3	27.9	18.2	58.0	0.000
- Aphasia subscale *	28.2	1.6	22.7	6.3	17.4	0.000
- Agnosia subscale *	24.9	2.7	16.8	8.2	33.5	0.000
- Dysexecution subscale	11.4	3.3	6.6	3.8	1.0	0.327
- Apraxia subscale	12.5	1.9	9.7	3.2	3.9	0.051
- Attention/ calculation subscale	10.1	4.8	4.7	4.8	0.3	0.612

* Mean scores of GDS stage 3 group are significantly different from the GDS stage 4 group ($p < 0.001$, by Levene's Test for Equality of Variances).

GDS: Global Deterioration Scale; CARDS-cog, Cognitive Assessment & Reference Diagnoses System-cognitive subscale; S.D., standard deviation

(adopted from Suh (2002), *Journal of the Korean Geriatric Psychiatry*)

General instruction for interviewer

To scrutinize this chapter would be a shortcut to be an excellent interviewer of the CARDS.

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How to deal with problems related to ageing and physical illnesses?

There's no need to consider one's age and/or presence of physical illnesses when you rate each item of the CARDS.

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Contents

The CARDS comprises a diagnostic test divided into the following sections. The psychometric test items in CARDS-cog were selected based on their sensitivity, brevity and ease of administration. The tests assessed verbal memory, concentration, language, perception, praxis and executive functioning, mainly focusing on the DSM-IV criteria for dementia diagnosis. The test is divided into a memory test and a non-memory test and is composed of 94 test/question/objective rating/full-up items.

● Neuropsychological Battery CARDS-cog

A. Memory test comprised 5 sections.

1. Word list memory – immediate recall: a 10-word list was presented and free recall was assessed across three trials, with the words being presented in a different order for each trial.
2. Word list memory – delayed recall: delayed memory for the 10 words in the word list (immediate recall) test was assessed approximately 10 minutes after the third immediate recall trial.
3. Word list memory –recognition: recognition of the 10 words from the word list (immediate recall) test was tested when the words were presented together with 10 distracter words.
4. Object list memory – delayed recall: delayed memory for the 10 objects in the object naming test was assessed approximately 10 min after the naming objects trial.

5. Object list memory – recognition: recognition of the 10 objects from the object naming test was tested when each was presented together with three distracter objects.
6. Recent memory: this was assessed by asking the names of the current and past Korean presidents and the food taken at the meal prior to examination. Questions about current news stories (within 1 or 2 weeks) and questions about their visit to the clinic (e.g. place he/she left for the clinic, time to reach the clinic, accompanying persons, transportation) were also asked.
7. Remote memory: this was assessed by asking about persons or events between 1910 and 1970 that would be very familiar to the Korean elderly. Two sets of tests were prepared to ask a different series of questions for subjects born before and after 1940.
8. Orientation: this was assessed by questioning the date, month, year, week day, season, name of the building, floor, area, city and person.
9. Today's weather: subjects were asked about the weather.
10. Current time: subjects were asked the current time.

B. Non-memory test comprised 5 sections.

1 Agnosia test:

- (a) object naming test : 10 3-Dimensional pictures of culturally familiar Korean objects were displayed and rotated so that they were shown on every side for subjects to name;
- (b) identification test : pictures are presented on a touch screen; (i) value of money: 10,000 won bill, 1,000 won bill and 500 won coin were presented for subjects to name the value; (ii) famous historical persons (King Sejong and Admiral Sushin Lee of Yi dynasty and the patriot Kwansoon Rhue in modern time) were presented for identification; (iii) essential signs for daily living: red traffic light, the sign for 'smoking prohibition' and the sign for 'toilet' were presented for identification; (iv) clock face reading: subjects were asked 'What time is it on the presenting clock?';
- (c) object constancy test : six pictures of objects taken from unusual angles were presented and subjects were asked 'What's this?'.

2 Aphasia test:

- (a) object naming test: this score was counted twice, both in the agnosia test and the aphasia test;
- (b) verbal fluency (animal category): this test assessed language and, more importantly, executive functioning by asking the subject to name as many animals with 4 legs as possible within 1 min;
- (c) word definition: this test assessed expressive language and semantic memory by asking the meaning of objects of two concrete nouns and an abstract noun; (
- (d) interviewer's judgment from free conversation; (i) expressive language: free production of speech; (ii) comprehension of spoken language: understanding during conversation; (iii) paraphasia and word-finding difficulty: pause or stammering.

3 Apraxia test:

- (a) ideational praxis;
- (b) ideomotor praxis;
- (c) constructional praxis: visual–spatial skills were assessed by presenting four line drawings of usual figures (a circle with wave inside, star, interlocking pentagons and flower in a pot) for subjects to copy;
- (d) subjects were asked to write a sentence.

4 Dysexecution test:

- (a) abstract thinking: differentiation and categorization of similar words;
- (b) ideational fluency: asking 'Please tell me how you can use water?';
- (c) visual reasoning: asking grouping strategies during oriental or western card games;
- (d) judgment: asking 'what you have to do in this situation?' and showing a drawing of a near-drowning man in the sea.

5 Attention/calculation test:

- (a) trail making test A: connect 15 points, forming a house shape and check reaction time;
- (b) serial 7 subtraction;
- (c) backward citation of weekdays;
- (d) combined arithmetic test in daily life (add/subtract /multiply/divide): usual examination at the level of 1st or 2nd grade in primary school.

- ***Interview with subject***

This section is composed of 15 questions related to depression, hallucination and delusion. Although most depressive symptoms are checked in this section, the rater's judgment is also considered in the diagnosis of depression. These questions were constructed to meet the DSM-IV criteria for depression.

- ***Interview with informants***

This section comprised 47 questions for the informant, a relative or a carer who knows the subject. It consisted of questions relating to cognitive difficulty, functional decline, activity of daily living, delirium, psychotic symptoms (delusions or hallucinations), time passed since the onset of symptoms, temporal relations between events (i.e. cerebrovascular accident and dementia). All questions were essential for the diagnosis of common mental disorders in later life and differential diagnosis for different diagnostic criteria. The questions were constructed with reference to a review of the clinical records of subjects with dementia to meet diagnostic criteria.

- ***Interviewer's judgment***

This section comprised five items for objective rating (i.e. depressed mood, flat affect, psychomotor retardation, perseveration, incoherent speech and behavior).

- ***clinician's examination***

This section comprised physical, neurological and mental state examination, including blood pressure, deep tendon reflexes, gait, extrapyramidal symptoms and paralysis of extremities or trunk.

- ***laboratory tests and brain imaging***

The results of laboratory, radiological and other investigations were recorded. Results were obtained from a complete blood count (CBC), liver function test (LFT) with urea and electrolyte, thyroid function test (TFT), vitamin B12 and folate, VDRL, EEG, CT/MRI and SPECT/PET.

PRACTICAL GUIDELINE FOR THE CARDS

The Developer of the CARDS, Dr. Suh, Guk-Hee is testing an elderly man with subjective memory impairment.

Neuropsychological test battery: CARDS-cog

CARDS-cog is a cognitive examination to test neuropsychological function by valid and reliable test battery of essential cognitive domains.

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The screenshot shows a software interface for a card-based memory test. On the left, there is a vertical column of five blue buttons: 'SKIP', 'NEXT', 'STOP', 'PREV', and 'QUIT'. The main area is divided into two sections: 'Question' and 'Answer'. The 'Question' section contains the text: '오늘 점심 식사를 하셨습니까?' and '식사하실 때 반찬이 무엇이었습니다?'. The 'Answer' section contains two radio buttons: '틀렸습니다' (Incorrect) and '맞았습니다' (Correct). At the bottom, there is a small text box with Korean instructions: '아침식사를 하지 않았다고 대답하면, 직전에 식사한 때를 기준으로 '어제 저녁 식사를 하셨습니까?' 라고 고쳐서 질문한다. 검사 현장에 있는 어느 누구도 객관적인 정답을 알 수 없는 상황이라면, 검사자 화면 왼쪽에 있는 'SKIP' BUTTON을 눌러야 한다. 대신 'NEXT' BUTTON을 눌러 다음 문제로 진행하면, 틀린 것으로 간주되어 부정확한 진단의 원인이 된다'.

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Serial 7 subtraction & Backward Citation of Weekdays

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The screenshot shows a software interface titled 'ED CARDS'. On the left side, there is a vertical column of five blue buttons: 'SKIP', 'NEXT', 'STOP', 'PREV', and 'QUIT'. The main area is divided into two sections: 'Question' and 'Answer'. The 'Question' section contains the text '일 요일 부터 시작하여 요일을 거꾸로 말해보세요.' (Start from the first day of the week and say the days of the week in reverse order). The 'Answer' section contains a row of seven checkboxes labeled '일', '토', '금', '목', '수', '화', and '월' (Sun, Sat, Fri, Thu, Wed, Tue, Mon). At the bottom of the interface, there is a small text box that reads: '앞에서 요일을 순서에 맞지 않게 말했지만 다음에 말한 요일부터는 순서에 맞으면 그 부분은 맞는 것으로 간주한다' (Even if you said the days of the week in the wrong order before, from the next day you say, if it matches the order, that part is considered correct).

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The screenshot shows the CARDs software interface. On the left side, there is a vertical column of buttons: SKIP, NEXT, STOP, PREV, and QUIT. The main area is divided into two sections: 'Question' and 'Answer'. The 'Question' section contains the text: '사과와 배의 공통점은 무엇입니까?'. The 'Answer' section contains a list of options with radio buttons: '틀렸습니다(예:동글다)', '음식,껍질,재배,영양가', and '과일'. At the bottom of the interface, there is a small text box with the following text: '면담 대상자가 ‘비슷한 점이 없다’고 말하면, ‘어딘가 공통점이 있습니다. 어떤 점이 비슷한지 말해주시겠습니까?’라고 말한다. 이 문항에 한해서만, ‘사과와 배의 공통점은 둘 다 과일이라는 점입니다.’라고 말해줄 수 있다'.

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
Laboratory tests and Brain Imaging

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Start

- To test New Clinet
- To restart at the item where you stopped

Report

- **Pentagonal Diagram and Diagnoses**
Crossing points of the red pentagon(test scores) nearer to center than those of the blue pentagon (cutoff points) mean dysfunctional among 6 key cognitive domains.
The computer-produced DSM-IV diagnosis is given in addition to the Reference diagnoses presenting all the possible diagnoses of the client to support your clinical diagnosis. In addition, MMSE score and modified Hachinski score are also given.
- **Summary of the CARDS-cog & Symptom Report**
Percentile profiles of the cognitive level and the gross features of the key symptoms are given.
- **Clinical, Laboratory & Radiological Report**
Key information of medical disease, summary of physical & neurological examination, and summary of the laboratory & radiological reports are given.
- **New & Old in the horizontal menu bar**
When you want to compare several CARDS trials in succession, select ‘New & Old’ in the horizontal menu bar.


Exercise

You are advised to select ‘Exercise’ to make a client familiar with computer and CARDS

Exit

Bye!

CHAPTER 3. PRACTICAL GUIDELINE FOR THE CARDS

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
Summary of the CARDS-cog & Symptom Report 이고

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THE PENTAGONAL DIAGRAM AND DIAGNOSES

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

SUMMARY OF THE CARDS-COG & SYMPTOM REPORT

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Symptom summary		
CATEGORY	CONTENTS	ANSWER
Amnesia	Forget where he/she put things	Severe
	Forget names	None
	Forget recent events	Severe
	Onset of Age	57.3 (year-old)
	Duration(unit:week)	52 (weeks)
Aphasia	Pattern of onset	Abrupt
Aphasia	Lose track while talking	Severe
Agnosia	Lost in a familiar place	Mild
	Lost in the house	None
Apraxia	Poor performance in complex task	Severe
Functional loss	IADL	37% (3 / 8)
	ADL	75% (6 / 8)
	Duration(unit:week)	52 (weeks)
	Pattern of functional loss	Abrupt or Stepwise
Disinhibition	Demoralizing, asocial behavior	Absent
Overeating	Especially sweet foods	Absent
Stereotypy	Repetitive behavior	Absent
	Repetitive speech	Absent
	Various features of stroke	Absent

CLINICAL, LABORATORY & RADIOLOGICAL REPORTS

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Past medical history	
CATEGORY	ANSWER
Hypertension	Absent
Myocardial Infarction	Absent
Parkinsonism	Absent
Head trauma	Absent
Seizure	Absent
Alcohol problem	Absent
Neuroleptic Malignant Syndrome	Absent
Physical and Neurological examination	
CATEGORY	ANSWER
Systolic Blood pressure	120 mmHg
Diastolic Blood pressure	080 mmHg
Deep tendon reflex	Normal
Babinski's reflex	Absent
Extrapyramidal symptom	Absent
Paralysis	None
Gait	Normal

Laboratory tests and brain imaging	
CATEGORY	Records
EKG	You may not order tests or enter results.
BUN/Cr,e'	You may not order tests or enter results.
LFT	You may not order tests or enter results.
Lipid battery	You may not order tests or enter results.
Vitamin B12	You may not order tests or enter results.
Folate	You may not order tests or enter results.
TFT	You may not order tests or enter results.
VDRL	You may not order tests or enter results.
HIV	You may not order tests or enter results.
Brain imaging	You may not order tests or enter results.

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A

abstract thinking, 13
ADAS-K-cog, 5
Agnosia test, iii, 4, 12
Aphasia test, iii, 4, 13
Apraxia test, iii, 4, 13
Attention/calculation test, 13
Attention/Calculation test, iii

B

backward citation of weekdays, 14

C

CAMCOG-K, 5
CARDS, iii, 82
CARDS-cog, iii, 2, 3, 4, 5, 6, 11, 17, 48, 59, 68, 75, 77
CBC, 15
cerebrovascular accident, 14
clinician's examination, iii, 5, 15
combined arithmetic test, 14
comprehension of spoken language, 13
constructional praxis, 13
CT/MRI, 15, 82
Current time, 12

D

deep tendon reflexes, 15
DSM-IV diagnosis, 5
Dysexecution test, iii, 4, 13

E

EEG, 15
Exercise, 8, 68
expressive language, 13
extrapyramidal symptoms, 15

F

F5 Key, 28
flat affect, 15, 58
folate, 15, 76

G

gait, 15, 67
GDS, 3, 6

I

ideational fluency, 13
ideational praxis, 13
identification test, 12
ideomotor praxis, 13
Interview with informants, iii, 4, 14, 59
Interview with subject, iii, 14, 48
interviewer's judgment, iii, 13
Interviewer's judgment, 4, 15, 54

J

judgment, 9, 13, 14

L

laboratory tests and brain imaging, iii, 5, 15
LFT, 15

M

Memory test, iii, 4, 11
MMSE, 1, 5, 68, 76

N

NEXT, 8
Non-memory test, iii, 12

O

object constancy test, 12
Object list memory – delayed recall, 11
Object list memory – recognition, 12
object naming test, 11, 12, 13
Orientation, 4, 12, 17

P

paralysis, 15
paraphasia, 13
perseveration, 15, 40
psychomotor retardation, 15

Q

Qualifying Interviewer, 10

R

Recent memory, 12
Reference diagnosis, 5
Remote memory, 12

S

serial 7 subtraction, 14
Skip, 7, 8
SPECT/PET, 15, 67

T

TFT, 15
Today's weather, 12
trail making test A, 14

V

VDRL, 15
verbal fluency, 13
visual reasoning, 13
vitamin B12, 15, 67

W

word definition, 13
Word list memory – delayed recall, 11
Word list memory – immediate recall, 11
Word list memory – recognition, 11
word-finding difficulty, 13
write a sentence, 13